


# C-Chip 4ch

CE

***Count your cells,  
for the most affordable price***



In the USA

MAYFLOWER  
BIOSCIENCE



[www.mayflowerbio.com](http://www.mayflowerbio.com)

Inquiries: [sales@mayflowerbio.com](mailto:sales@mayflowerbio.com)

Support: [support@mayflowerbio.com](mailto:support@mayflowerbio.com)

Phone: 314-485-5210




# C-Chip 4ch

- ✓ New 4-channel type disposable hemocytometer
- ✓ Created by the best-known microfluidic slide manufacturer
- ✓ High Quality, Low Price

## C-Chip 4channel is a Disposable Plastic Hemocytometer

- No need for coverslips
- Disposable, eliminates the need for washing
- Reduces exposure to infectious samples and risk of contamination
- Guaranteed reproducibility and reliability
- Four tests
- Made of sturdy quartz grade optical plastic




## Models

No.	Grid	Depth
DHC-N04	Neubauer Improved	100 $\mu$ m

## Applications

Blood analysis (Hematology) : Blood cell counting  
 Cell Culture : Cell concentration measurement / Cell viability test  
 Microbiology : Bacteria & Fungal spore counting  
 IVF,IUI : Sperm counting



## Specifications

Dimension : 25 (W) x 75 (L) x Thickness 1.6 mm  
 Chamber depth : 100  $\mu$ m  
 Chamber volume : 10  $\mu$ L

In the USA

**MAYFLOWER**  
 BIOSCIENCE



[www.mayflowerbio.com](http://www.mayflowerbio.com)

Inquiries: [sales@mayflowerbio.com](mailto:sales@mayflowerbio.com)

Support: [support@mayflowerbio.com](mailto:support@mayflowerbio.com)

Phone: 314-485-5210